


VITROCELL® in vitro Exposure Technology

Aerosol exposure devices


for *in vitro* toxicity testing at the air/liquid interface




Aerosol exposure devices for *in vitro* toxicity testing at the air/liquid interface

Mission: use of turnkey *in vitro* exposure systems incl. online dosimetry tools

Gases Complex Mixtures Airborne Particles (Nano Particles)




source: environment.scotland.gov.uk



source: cysticfibrosisnewstoday.com

Cell Types of Respiratory Tract

Function	Zone	respiratory tract scheme	Anatomy	Airway Surface	No. Of Airways	Air space/Liquid layer	Particle Diameter	Particle size	Cell Types
Air conditioning, temperature and humidity, cleaning. Fast particle clearance, air conduction	Extrapulmonary	 <small>source: W. Koch, Fraunhofer ITEM, Hannover</small>	Nose, Mouth, Esophargus	$2.4 \times 10^{-3} \text{m}^2$	1	18 mm/ 10 μm	5-10 μm	Dust, Mist, fly-ash, Pollen, Spores	Ciliated Cells, non-ciliated cells (Globet/Mucous/Serous/Brush/Endocrine/Basal/Intermediate Cells)
	Pulmonary		Trachea	$4.5 \times 10^{-2} \text{m}^2$	1		3-5 μm	Fog, Mold, Bacteria, Exhaust gases	
			Main Bronchi/Bronchi	$3 \times 10^{-2} \text{m}^2$	511		2-3 μm		
Air conditioning, gas exchange, slow particle clearance			Bronchioles/terminal bronchioles	$2.6 \times 10^{-4} \text{m}^2$	$6,5 \times 10^4$	0.5-1mm/ 3 μm	1-2 μm	Respiratory epithelium (ciliated and non-ciliated cells with clara cells (without globet cells)	
			Respiratory bronchioles	7.5m^2	$4,6 \times 10^5$			Respiratory epithelium with maily secretory cells and few ciliated cells	
gas exchange, very slow particle clearance				Alveolar ducts and sacs	140m^2	$4,5 \times 10^7$	0.4mm/ 0.1 μm	0.01-1 μm	Fumes, Tobacco smoke, viruses, Oil smoke

Source: Encyclopaedia of Occupational Health and Safety , Chapter 10 - Respiratory System



Costar® Thincert® Falcon®

Examples of cell culture inserts

Human Cell Culture Models

Airway epithelial cells

Calu-3 (ATCC HTB-55)
16HBE14o- (can be obtained from D.C. Gruenert)
BEAS-2B (ATCC CRL-9609)
NuLi-1 (ATCC CRL-4011)

Primary airway epithelial cells

hBEpC

Commercially available from Epithelix (CH), MatTek (USA)

Alveolar epithelial cell lines

A549 (ATCC CL-185): ATII phenotype
Immortalized human ATII cells with ATI phenotype
NCI-H441 (ATCC HTB-174): ATII and Clara cell phenotype

Commercially available from inscreenex (D)

Primary alveolar epithelial cells

hAEP: ATII cells that differentiate in vitro into ATI-like morphology

3D cultures

3D aggregates of A549 cells
Bilayer co-culture model: epithelial & endothelial cells
Triple cell co-culture model: epithelial cells, macrophages, dendritic cells
Double, triple and quadruple cell co-culture models: epithelial cells, endothelial cells, mast cells, macrophages



Costar® Thincert® Falcon®

Examples of cell culture inserts

AT = alveolartype

Source: Barbara Rothen-Rutishauser, Advanced in vitro lung models in nanotoxicology research – advantages and limitations, Adolphe Merkle Institute
University of Fribourg Marly, Switzerland
Jud et al. Swiss Med Wkly 143:0 (2013).

Human Cell Culture Models

Airway epithelial cells

Calu-3 (ATCC HTB-55)
16HBE14o- (can be obtained from D.C. Gruenert)
BEAS-2B (ATCC CRL-9609)
NuLi-1 (ATCC CRL-4011)

Primary airway epithelial cells

hBEpC

Alveolar epithelial cell lines

A549 (ATCC CL-185): ATII phenotype
Immortalized human ATII cells with ATI phenotype
NCI-H441 (ATCC HTB-174): ATII and Clara cell phenotype

Primary alveolar epithelial cells

hAEpC: ATII cells that differentiate in vitro into ATI-like morphology

3D cultures

3D aggregates of A549 cells
Bilayer co-culture model: epithelial & endothelial cells
Triple cell co-culture model: epithelial cells, macrophages, dendritic cells
Double, triple and quadruple cell co-culture models: epithelial cells, endothelial cells, mast cells, macrophages

Biomimetic microsystems

Breathing lung-on-a-chip: epithelial & endothelial cells
Perfused chip: epithelial & endothelial cells

AT = alveolartype

Source: Barbara Rothen-Rutishauser, Advanced in vitro lung models in nanotoxicology research – advantages and limitations, Adolphe Merkle Institute
University of Fribourg Marly, Switzerland
Jud et al. Swiss Med Wkly 143:0 (2013).

Commercially available from Epithelix (CH), MatTek (USA)

Commercially available from inscreenex (D)

Alveolix (CH), Emulate (USA)



Costar® Thincert® Falcon®

Examples of cell culture inserts




AlveoliX
Cooperation VITROCELL




Emulate

Disadvantages Submerged Exposure



Submerged Exposure
Low sensitivity due to contact
with media
> Undefined Dose




➡ Different sedimentation times


➡ Upward diffusion

➡ Buoyant behaviour


Advantages of Air/Liquid Exposure



Low Sensitivity
Undefined Dose





Low Sensitivity
Undefined Dose



High Sensitivity
Defined Dose
Physiologically Relevant

How to test airborne substances *in vitro* at the air/liquid interface?



Steps for ALI Exposure

- Provide the Test Substance

- Aerosol generation (Lab)
- Direct sampling (Environment)

Substance Characterisation

- Define Exposure Principle & Device

- > Continuous Flow
- > Single Droplet Sedimentation
- > Dry Powder

- Establish Dose / Response Range

- Dosimetry Tools

- Gas phase
- Particle phase

Critical System Elements for ALI Exposure

- Provide the Test Substance

- Aerosol generation (Lab)
- Direct sampling (Environment)

Substance Characterisation

- Define Exposure Principle & Device

- > Continuous Flow
- > Single Droplet Sedimentation
- > Dry Powder

- Establish Dose / Response Range

- Dosimetry Tools


- Gas phase
- Particle phase

Providing the Test Substance > in the Lab

Gases & Complex Mixtures



Gas Cylinder




Smoking Robot



Automotive Engine

Liquid Aerosols


Dry Powders



Aerosol Generator



Cloud Aerosol Generation



Bulk Powder Aerosol Generation



Small Quantity Powder Chamber

Providing the Test Substance > Sampling directly from Environment




Gases & Complex Mixtures



Mobile Test Lab



Automated Exposure Station
> Suitable for long-term exposures (24 h)



Liquid Aerosols

Dry Powders

Critical System Elements for ALI Exposure

- Provide the Test Substance

- Aerosol generation (Lab)
- Direct sampling (Environment)

Substance Characterisation

- Define Exposure Principle & Device


- > Continuous Flow
- > Single Droplet Sedimentation
- > Dry Powder

- Establish Dose / Response Range

- Dosimetry Tools

- Gas phase
- Particle phase

Exposure Principle 1: Continuous Flow




- Aerosol maintains reactivity of original compounds
- Physiologically relevant
- Possibility to integrate online dosimetry tools

Air/Liquid Interface
High Sensitivity
Defined Dose


Continuous Flow Device - VITROCELL® 12/6 CF

VITROCELL® 12/6 CF exposure module rack for 12- and 24- well sized inserts



Continuous Flow Device - VITROCELL® 6/4 and 6/3


VITROCELL® 6/4 CF and 6/3 CF exposure modules for 6-, 12- and 24- well sized inserts





Optional Standalone Deposition Enhancement



- Working by electrical field
- Approx. 6....8-fold increase of deposition



Number size distribution of fluorescein sodium in the ALI aerosol measured by SMPS (mean of 29), $d_{\text{modal}} = 65 \text{ nm}$, $s_g = 2.0$



Data courtesy of KIT, Karlsruhe, Germany


Automated Exposure Stations




- Automated process control
- System is entirely heated (37°C)
- Integrated humidification of aerosol and clean air (85% r.H.)
- Aerosol flow control by mass flow controllers
- Online dose measurement (QCM)

Exposure Principle 2: Cloud Single Droplet Sedimentation

- Exposure systems for small quantities of liquid aerosols (200 µl)
- High droplet output rate – Cloud dynamics with high deposition efficiency
- Easy-to-use and clinically relevant




Exposure Principle 2: Cloud Single Droplet Sedimentation



Exposure Principle 2: Cloud Single Droplet Sedimentation


- Exposure systems for small quantities of liquid aerosols (200 µl)
- High droplet output rate – Cloud dynamics with high deposition efficiency
- Easy-to-use and clinically relevant




Versions for 6-, 12- and 24-well sized inserts



Well-to-well variability between inserts
Range 3 and 10% (1 σ)




Exposure Principle 3: Small Quantities of Dry Powder / Powder Chamber



Powder Chamber Process

Sedimentation Tubes

- Variable length (regular: 200 mm)
- Optional heating for Thermophoresis, to decrease sedimentation time



Exposure Tray

- Heated to 37 °C

Step 1: Filling and Sedimentation Time – usual mass of 20 mg per exposure

Substances / Preferred Exposure Principles

Substances



Gases



Environmental Atmospheres



Vehicle Emissions



Cigarette Smoke & e-Cig Vapors



Particles & Nanomaterials



Medical Liquid Aerosols



Powders

Exposure Principle / Module Type

Continuous Flow



Continuous Flow



Continuous Flow



Continuous Flow



Continuous Flow



Cloud



Powder Chamber



Cloud



Continuous Flow



Powder Chamber



Critical System Elements for ALI Exposure

- Provide the Test Substance

- Aerosol generation (Lab)
- Direct sampling (Environment)

Substance Characterisation

- Define Exposure Principle & Device

- > Continuous Flow
- > Single Droplet Sedimentation
- > Dry Powder



- Establish Dose / Response Range

- Dosimetry Tools

- Gas phase
- Particle phase

Dose/Response Relations

- Time-based exposure at same concentration
- High to low exposure concentration levels



Critical System Elements for ALI Exposure

- Provide the Test Substance

- Aerosol generation (Lab)
- Direct sampling (Environment)

Substance Characterisation

- Define Exposure Principle & Device


- > Continuous Flow
- > Single Droplet Sedimentation
- > Dry Powder

- Establish Dose / Response Range

- Dosimetry Tools

- Gas phase
- Particle phase

System Element: Online Dosimetry Tools




Online dosimetry hotspots



Dose Monitoring with Microbalance Technology



Please visit our website vitrocell.com and the „Publications“ chapter with over 170 documents



Advanced in vitro exposure systems.

[INHALATION TOXICOLOGY](#)[SKIN EXPOSURE](#)[AUXILIARY EQUIPMENT](#)[INFORMATION CENTER](#)

Working together?
Partnership for success!

Valuable documentation on all VITROCELL® products, references as well as a compilation of relevant publications make this section a powerful information source for in vitro research.

[Downloads](#)[News](#)[Publications](#)[References](#)[Suppliers](#)[About us](#)[Useful Links](#)[Contact](#)

[Home](#) | [News](#)

10. May. 2018

Combining MucilAir™ and Vitrocell® Powder Chamber for the In Vitro Evaluation of Nasal Ointments in the Context of Aerosolized Pollen

<https://doi.org/10.3390/pharmaceutics10020056>

Julia Metz ^{1,2}, Katharina Knoth ¹, Henrik Groß ¹, Claus-Michael Lehr ^{1,2,3}, Carolin Stäbler ⁴, Udo Bock ⁵, and Marius Hittinger ¹

¹Department of Drug Delivery, PharmBioTec GmbH, 66123 Saarbrücken, Germany





²Department of Biopharmaceutics and Pharmaceutical Technology, Department of Pharmacy, Saarland University, 66123 Saarbrücken, Germany

³Helmholtz Institute for Pharmaceutical Research Saarland (HIPS), Helmholtz Center for Infection Research (HZI), 66123 Saarbrücken, Germany

⁴Department of Scientific Affairs Consumer Health, Bayer Vital GmbH, 51368 Leverkusen, Germany

⁵Bock Project Management, 54456 Tawern, Germany

CATEGORIES

Exhibitions (33)	
Product News (41)	
Publications (171)	
User Group Meeting (8)	

TAGS

aerosol generator aerosolist
airborne particles
air-liquid Interface AMES
autosampler congress diesel
e-cigarettes exposure station

Thank you for your attention !

